

MIDLABS®



Providing Quality & Innovative
**Ophthalmic Surgical
Systems & Instruments
Since 1982**



www.midlabs.com

Made in the U.S.A. • ISO 13485:2016 Certified

MID Labs: Elevating Vision Care

Nestled in the vibrant San Francisco Bay Area, MID Labs stands as a beacon of innovation in ophthalmic healthcare. Our resolute mission is to exceed expectations by delivering ingenious, world-class solutions that safeguard and enhance the vision of patients across the globe.

Our exceptional team of dedicated professionals embodies the values of collaboration, initiative, integrity, and uncompromising attention to precision.

At MID Labs, we are more than innovators; we're visionaries. As pioneers of vitrectomy technology, we have revolutionized surgical procedures, crafting instruments that set global benchmarks. Our diverse product line, distributed in over 30 countries and continuously expanding, features Vitreous Cutters, Bi-Blade Vitreous Cutters, and the Single-Step Cannula Insertion System.

Select MID Labs as your trusted partner.

**Together, we are not merely
shaping the future of eye health;
we are illuminating it.**

Our Business Throughout the Years

- 1981** MID Labs is founded by Dr. Carl Wang. Dr. Wang was associated with vitreoretinal surgery from its infancy. In 1975, his Ocutome® 800 became the first complete system for posterior vitrectomy.
- 1982** The MicroVit® system introduced venturi aspiration to vitreoretinal surgery and the MicroVit® probe provided the first disposable high-quality vitreous cutter surgical system.
- 1994
1997** MID Labs introduced the SupraVit® VitreoRetinal Surgical System with a revolutionary aspiration system and an Extended Use Vitreous Cutter.
- 2000
2002** MID Labs launched the Vit Enhancer® and Adaptable Vit Enhancer® as High Speed Vitreous Cutting Systems, creating an evolution in vitreous cutting.
- 2005** MID Labs introduced a non-system exclusive 25 Gauge Vitreous Cutter line, making micro-incision vitrectomy available for a wide variety of vitreoretinal surgical systems.
- 2009** MID Labs introduced the 27 Gauge Vitreous Cutters, 25 and 23 Gauge Cannula Insertion Systems as well as developing specialty Vitreous Cutters, for OEM partners.
- 2011** MID Labs introduces the UVE (Ultimate Vit Enhancer®) offering an 8000 cpm capability, allowing enhanced vitreous cutting capabilities to be added to existing surgical systems.
- 2014** MID Labs receives Joint Venture financing led by OrbiMed and Frontline BioVentures and is well positioned to leverage its technology and leadership in vitrectors while broadening its product lines on a global basis.
- 2016** MID Labs introduces Bi-Blade®. The new MID Labs Bi-Blade® Vitrectomy Cutter provides enhanced performance, improved surgical control, and increased surgical efficiency.
- TODAY** MID Labs continues to innovate with its own product line and is a supplier to major ophthalmic system manufacturers in Europe, Asia and within the USA.

MID Labs Vision

"To Excel at providing Quality and Innovative Ophthalmic Products that Benefit the Patient"

CONTACT US

When placing an order please refer to the catalog number on the backside of this brochure.

✉ sales@midlabs.com

☎ +510-357-3952

🇺🇸 1-800-929-5227

📞 +510-357-1582

WHERE THERE IS VISION, THERE IS HOPE

557 McCormick St. San Leandro, CA 94577
www.midlabs.com

CE0197




LM00013-Rev A

MID Labs Products


Bi-Blade® and Conventional Vitreous Cutters





Below products are sold 6 per box, sterile with tubing sets

Photo and Description	Catalog Number
 Vitreous Cutter: 13.2 psi, BP Extended PVC Tubing Set, Sterile, Single Use	• 2500CE = (20 gauge) • 2520CE = (23 gauge) • 2540CE = (25 gauge)

Below products are sold 6 per box, sterile with tubing sets



Photo and Description	Catalog Number
 Vitreous Cutter: 25 psi, Luer Posterior PVC Tubing Set, Sterile, Single Use	• 2400CE = (20 gauge) • 2420CE = (23 gauge) • 2440CE = (25 gauge)

Below products are sold 6 per box, sterile with tubing sets


Photo and Description	Catalog Number
 Vitreous Cutter: 25 psi, Luer Anterior PVC Tubing Set, Sterile, Single Use	• 2450CE = (20 gauge) • 2470CE = (23 gauge) • 2490CE = (25 gauge)
 BB Vitreous Cutter: 25 psi, Luer Anterior PVC Tubing Set, Sterile, Single Use	• 2470QB = (23 gauge) • 2490QB = (25 gauge)

MID Labs Products



Below products are sold 6 per box, sterile with tubing sets

Photo and Description	Catalog Number
 Vitreous Cutter: 25 psi, Luer Extended PVC Tubing Set, Sterile, Single Use	• 2700CE = (20 gauge) • 2720CE = (23 gauge) • 2740CE = (25 gauge)
 BB Vitreous Cutter: 25 psi, Luer Extended PVC Tubing Set, Sterile, Single Use	• 2720BB = (23 gauge) • 2740BB = (25 gauge) • 2760BB = (27 gauge)


Below products are sold 6 per box, sterile with tubing sets

Photo and Description	Catalog Number
 Vitreous Cutter: 27.5 psi, CPC Posterior PVC Tubing Set, Sterile, Single Use	• 2900CE = (20 gauge) • 2920CE = (23 gauge) • 2940CE = (25 gauge)


Below products are sold 6 per box, sterile with tubing sets

Photo and Description	Catalog Number
 Vitreous Cutter: 20 psi, Luer Anterior PVC Tubing Set, Sterile, Single Use	• 2350CE = (20 gauge) • 2370CE = (23 gauge) • 2390CE = (25 gauge)
 BB Vitreous Cutter: 20 psi, Luer Anterior PVC Tubing Set, Sterile, Single Use	• 2370BB = (23 gauge) • 2390BB = (25 gauge) • 2310BB = (27 gauge)

Below products are sold 6 per box, sterile with tubing sets

Photo and Description	Catalog Number
 BB Vitreous Cutter: 13.2 psi, Luer Posterior PVC Tubing Set, Sterile, Single Use	• 2320BB = (23 gauge) • 2340BB = (25 gauge) • 2360BB = (27 gauge)


Below products are sold 6 per box, sterile with tubing sets

Photo and Description	Catalog Number
 BB Vitreous Cutter: 37.5 psi, CPC Posterior PVC Tubing Set, Sterile, Single Use	• 2820BB = (23 gauge) • 2840BB = (25 gauge) • 2860BB = (27 gauge)

MID Labs Products

Cannula Insertion Systems

- ✓ Single Step
- ✓ Measuring and Marking Feature
- ✓ High-Flow Infusion




Below products are sold 6 per box, sterile with tubing sets

Item	Description	Catalog Number
CIS 23 Gauge	Cannula Insertion System: 23 Gauge, Sterile, Single Step, Gen II with Tubing Set	• 5700GEN • 5700GST
	CIS 25 Gauge	Cannula Insertion System: 25 Gauge, Sterile, Single Step, Gen II with Tubing Set
CIS 27 Gauge	Cannula Insertion System: 27 Gauge, Sterile, Single Step, Gen II with Tubing Set	• 5500GEN • 5500GST

When placing an order please refer to the catalog number on the backside of this brochure.

- ✉ sales@midlabs.com 🇺🇸 1-800-929-5227
- ☎ +510-357-3952 📠 +510-357-1582

The products in this brochure may contain information or intended uses that are not approved or offered in your country. Product specifications are subject to modification in design and scope due to ongoing technological advancements. For more details, please contact our Sales team.



MID Labs Products Coming Soon

Anterior & Posterior Vitreous Cutter Assemblies


Photo and Description	Catalog Number
 BB Vitreous Cutter: 37.5 psi, CPC Anterior PVC Tubing Set, Sterile, Single Use	• 2870BB = (23 gauge) • 2890BB = (25 gauge) • 2810BB = (27 gauge)



Photo and Description	Catalog Number
 BB Vitreous Cutter: 27.5 psi, CPC Anterior PVC Tubing Set, Sterile, Single Use	• 2970BB = (23 gauge) • 2990BB = (25 gauge) • 2910BB = (27 gauge)

Photo and Description	Catalog Number
 BB Vitreous Cutter: 31 psi, CPC Posterior PVC Tubing Set, Sterile, Single Use,	• 2620BB = (23 gauge) • 2640BB = (25 gauge) • 2660BB = (27 gauge)